



Data Support Co., Inc.

Fat Tester Comparison Chart

800-726-5883

www.dsctest.com

Benefits of the Digital DSC HFT2000F & DSC50F over Traditional "Drip" Methods

Featured Models:

HFT2000F®

DSC50F®

Hobart F-101®

Univex FA73®



Accuracy	± 1% Fat	± 1% Fat	± 3% Fat	± 4% Fat
Fat Range Measure	1-99%	1-99%	10-40%	10-40%
Test Time	15 Minutes	15 Minutes	15 Minutes	9 Minutes (3 Test Required)
PC Interface	Yes	Yes	No	No
Printer Options	Yes	Yes	No	No
Consumables	Disposable (No Clean Up)	Disposable (No Clean Up)	Wash & Re-Use	Wash & Re-Use
Unit Gross Weight	9lbs	12lbs	21lbs	12lbs
Effect of Initial Sample	No Effect	No Effect	Effects Results	Volume vs. Weight (Very Inaccurate)

The DSC Line of Fat Testers allows control of fat content at profitable levels and is cost effective for both small meat markets and large Chain Super Stores. This system is used in hundreds of markets with proven reliability and a fast pay-back period.